



AccuFix™ Cantilever Intracranial, Head & Neck Device RT-4525L



Patented



Accessories:



Traditional, Slimline &
Split-Frame Masks
(see page 38 & 39 for more sizes
and descriptions)



RT-4525TILT
Tilt-Pro Tilting Base
(see page 40 for more information)



RT-4525LSPLIT
Split Frame Adapter

This patented head and neck device retracts the shoulders away from the neck. The device's unique cantilever design allows for 360-degree treatment with minimal attenuation. The fully adjustable carbon fiber Shoulder-Loc™ shoulder depression system is designed to increase setup repeatability and reduce patient motion without the need for large sheets of thermoplastic over the shoulders.

- Adapters are available to mount the AccuFix™ to treatment couch tops and position the device in a fixed and repeatable location to prevent tipping that can occur with other cantilevered head and neck devices.
- A Wall Mount is available for the AccuFix™, shown on page 122.



RT-4525KVL
kVue™ AccuFix™
Cantilever Insert



RT-4525LMRI
AccuFix™ Cantilever MR
Head & Neck Device



RT-4525L	AccuFix™ Cantilever Head & Neck Device
RT-4525KVL	kVue™ AccuFix™ Cantilever Insert
RT-4525LMRI	AccuFix™ Cantilever MR Head & Neck Device
RT-4525TILT	Tilt-Pro Tilting Base
RT-4525LSPLIT	Split Frame Adapter

Shim is included with order of RT-4525L, additional shims are available for order (RT-4529S)